



Product Name: Disposable High-pressure Syringes



STERILE EO



MRI Syringe Kits

60ml/60ml Syringe

*250cm T Connector*Spikes

P/N : 300106

For use with Nemoto CT Dual shot Injectors
Alpha 7 CT

Description:

- 2*60ml Syringes
- 1*2500mm(96") T Coiled Tube 350PSI
- 1* Large Spike 1* Small spike



Components / Materials

Luer Type

- Male Luer Connector
- Female Luer connector

Materials

- Polyethylene Terephthalate (PET)
- Synthetic Rubber
- Polycarbonate (PC)
- Polyvinyl Chloride (PVC)—TOTM
- / Polyvinyl Chloride (PVC)—DEHP
- Polyethylene (PE)

Exclusive Distributed by:

Shenzhen Medis Medical Supplies Co., Ltd.

Rm1001, Block A, Building 3, Tianan Cloud Park, Bantian Street, Longgang District, Shenzhen 518129, China.

Made in China

www.medismed.com

Package / Sterilization

Package

- Blister package—Each kit
- Paper box package—Outer box
- Quantity per case—50kits
- Quantity per tray —1250kits
- Lid — Tyvek Medical paper
- Sterilization — EO GAS
- Gross Weight per case: 10kg
- Case size : 57*39*39cm

Shelf Life

3 Years

For more distributor information,
please contact medismed direct
medis@medismed.com

MEDIS
For Healthy Life