



Product Name: Disposable High-pressure Syringes



## MR Syringe Kits

### 60ml Syringe \*250cm Connector \*spikes



STERILE EO



P/N : 500102

For use with Medtron MR Injectors

### Accutron MR

#### Description:

2\*60ml Syringes

1\*2500mm(60") T Coiled Tube 350PSI

1\*Large Spike 1\* Small Spike

#### Components / Materials

##### Luer Type

Male Luer Connector

Female Luer connector

##### Materials

Polyethylene Terephthalate (PET)

Synthetic Rubber

Polycarbonate (PC)

Polyvinyl Chloride (PVC)—TOTM

/ Polyvinyl Chloride (PVC)—DEHP

Polyethylene (PE)

Exclusive Distributed by:

**Shenzhen Medis Medical Supplies Co., Ltd.**

Rm1001, Block A, Building 3, Tianan Cloud Park, Bantian Street, Longgang District, Shenzhen 518129, China.

Made in China

[www.medismed.com](http://www.medismed.com)

#### Package / Sterilization

##### Package

Blister package—Each kit

Paper box package—Outer box

Quantity per case—50kits

Quantity per tray —1250kits

Lid — Tyvek Medical paper

Sterilization — EO GAS

Gross Weight per case: 10kg

Case size : 57\*39\*39cm

##### Shelf Life

3 Years

For more distributor information, please contact medismed direct

[medis@medismed.com](mailto:medis@medismed.com)