



ASA-200E



Shenzhen Anke High-tech Co., Ltd.
Address: Block B, LingYa Industrial Zone, Tangtou
No.1 Road,
Bao'an District, Shenzhen, 518108, P.R.China
Tel: +86-755-21622518 26688889
Fax: +86-755-26695307 26685908
Email: anke@anke.com Skype: anke.1986
Website: www. anke.com



ANKEIJ-2020-V.2.0-ENG



The new structure, lightweight and flexible, saving you time and effort; Well-designed syringes of both A and B capacities are not only suitable for adult CT enhancement, but also for children's CT enhancement and angiography;

Staged injections to ensure enhanced CT results; Rich input and output interfaces can be linked with various CT models; 150cm connecting tube can secure continuous injection during the scanning process;

Dual consoles respectively for scanning room and operating room

100 sets of default injection programs, easily modifiable at any time use

'Air exhaust' confirmation button ensures the safety of the injection
Injecting interruption function handles emergency situations flexibly



ASA-200E Specifications:

Barrel Capacity	200ml
Injection Volume	1—200ml
Flow Rate	0.1—8ml/s (increment 0.1ml/s)
Interval Time	0—240s
Delay Time	0—99s
Default Programs	100 sets
Pahses	1—4
Pressure Limit	100psi、150psi、200psi、250psi
Power	220VAC±10%/50Hz±1Hz

Signal output interface provided, can be used with various types of CT scanners

