





CT-MRI-Angio Contrast Media Syringes





Gallery

Dedicated

Open

Honest

Sincerely

Profession

Friendly


Fast

Reliable

Patient

Trustable


MEBIS
For Healthy Life



Overview




1. Medrad Syringes-----Page 1-2
2. LF Syringes -----Page 3
3. Nemoto Syringes -----Page 4-5
4. Bracco Syringes-----Page 6
5. Medtron Syringes -----Page 7
6. Imaxeon Syringes -----Page 8
7. Anke/Sino/Seacrown Syringes-----Page 9-10
8. Tubings & Spike& Transfer Sets & Ulrich Tubings-----Page 11-12




P/N : 100101
Medrad MCT & MCT plus CT,
Vistron CT, Envision CT.

1*200ml Syringe 1* 1500mm Coiled Tube 1* J Tube (QFT)




P/N : 100101S
Medrad MCT & MCT plus CT,
Vistron CT, Envision CT.

1*200ml Syringe 1* 1500mm Coiled Tube 1* Large Spike




P/N : 100102
Medrad Stellant CT Single

1*200ml Syringe 1* 1500mm Coiled Tube 1* Large Spike




P/N : 100102Q
Medrad Stellant CT Single

1*200ml Syringe 1* 1500mm Coiled Tube 1* J Tube (QFT)



P/N : 100102D
Medrad Stellant CT Dual

2*200ml Syringes 1* 1500mm Y Coiled Tube 2*Spikes



P/N : 100102QD
Medrad Stellant CT Dual

2*200ml Syringes 1* 1500mm Y Coiled Tube 2* J Tubes (QFTs)



P/N : 100103
Medrad Spectris MR

2*65ml Syringes 1* 2500mm T Coiled Tube
1* Small Spike 1*Large spike



P/N : 100104
Medrad Spectris Solaris MR

1*65ml Syringe 1*115ml Syringe 1* 2500mm T Coiled Tube
1* Small Spike 1*Large spike



P/N : 100104N
Medrad MRXPerion MR

1*65ml Syringe 1*115ml Syringe 1* 2500mm T Coiled Tube
1* Small Spike 1*Large spike



P/N : 100105
Medrad Mark V & Mark V provis Angio.

1*150ml Syringe 1* J Tube (QFT)



P/N : 100106
Medrad Mark V & Mark V provis Angio.

1*60ml Syringe 1* J Tube (QFT)



P/N : 100107
Medrad Mark III & Mark IV Angio.

1*130ml Syringe 1* J Tube (QFT)



P/N : 100108
Medrad Mark VII & Avanta Angio.

1*150ml Syringe 1* J Tube (QFT)





P/N : 200101
Liebel-Flarsheim CT9000 & CT9000 ADV
CT, Optivantage CT, Optione CT.

1*200ml Syringe 1*1500mm Coiled Tube 1* J Tube (QFT)



P/N : 200101S
Liebel-Flarsheim CT9000 & CT9000 ADV
CT, Optivantage CT, Optione CT.

1*200ml Syringe 1*1500mm Coiled Tube 1* Large Spike




P/N : 200101D
Liebel-Flarsheim CT9000 & CT9000 ADV
CT, Optivantage HDD CT.

2*200ml Syringes 1* 1500mm Y Coiled Tube 2* J Tubes (QFTs)




P/N : 200101SD
Liebel-Flarsheim CT9000 & CT9000 ADV
CT, Optivantage HDD CT.

2*200ml Syringes 1* 1500mm Y Coiled Tube 2*Spikes



P/N : 200102
Liebel-Flarsheim Optistar MR, Elite MR.

2*60ml Syringes 1* 2500mm T Coiled Tube
1* Small Spike 1*Large spike




P/N : 200103
Liebel-Flarsheim Angiomat Illumena
Neo Angio.

1*150ml Syringe 1* J Tube (QFT)




P/N : 200104
Liebel-Flarsheim Angiomat 6000 Angio.

1*150ml Syringe 1* J Tube (QFT)




P/N : 300101
Nemoto A25 & A60 & A300 , Smart shot Alpha, Dual Shot GX, Dual Shot Alpha CT.

1*200ml Syringe 1*1500mm Coiled Tube 1* J Tube (QFT)



P/N : 300101S
Nemoto A25 & A60 & A300 , Smart shot Alpha, Dual Shot GX, Dual Shot Alpha CT.

1*200ml Syringe 1*1500mm Coiled Tube 1* Large Spike




P/N : 300101D
Nemoto Dual Shot, B200 CT.

2*200ml Syringes 1* 1500mm Y Coiled Tube 2* J Tubes (QFTs)



P/N : 300101SD
Nemoto Dual Shot, B200 CT.

2*200ml Syringes 1* 1500mm Y Coiled Tube 2*Spikes



P/N : 300102
Nemoto A25 & A60 & A300 CT, Smart Shot Alpha CT, Dual Shot CT.

1*100ml Syringe 1*1500mm Coiled Tube 1* J Tube (QFT)



P/N : 300102S
Nemoto A25 & A60 & A300 CT, Smart Shot Alpha CT, Dual Shot CT.

1*100ml Syringe 1*1500mm Coiled Tube 1* Large Spike



P/N : 300102D
Nemoto Dual Shot CT.

2*100ml Syringes 1* 1500mm Y Coiled Tube 2* J Tubes (QFTs)



P/N : 300102SD
Nemoto Dual Shot CT.

2*100ml Syringes 1* 1500mm Y Coiled Tube 2*Spikes



P/N : 300103
Nemoto Dual Shot Alpha CT.

1*200ml Syringe 1*100ml Syringe
1* 1500mm Y Coiled Tube 2* J Tubes (QFTs)



P/N : 300103S
Nemoto Dual Shot Alpha CT.

1*200ml Syringe 1*100ml Syringe
1* 1500mm Y Coiled Tube 2*Spikes



P/N : 300104
Nemoto Dual Shot GX CT.

1*200ml Syringe 1*60ml Syringe
1* 1500mm Y Coiled Tube 2* J Tubes (QFTs)



P/N : 300104S
Nemoto Dual Shot GX CT.

1*200ml Syringe 1*60ml Syringe
1* 1500mm Y Coiled Tube 2*Spikes



P/N : 300105
Nemoto Dual Shot CT.

1*100ml Syringe 1*60ml Syringe
1* 1500mm Y Coiled Tube 2* J Tubes (QFTs)



P/N : 300105S
Nemoto Dual Shot CT.

1*100ml Syringe 1*60ml Syringe
1* 1500mm Y Coiled Tube 2*Spikes



P/N : 300106
Nemoto Sonic Shot GX & Shot 7,
& Shot 50 MR.

2*60ml Syringes 1* 2500mm T Coiled Tube
1* Small Spike 1*Large spike



P/N : 300107
Nemoto Nemoto 120s & Press Pro Angio.

1*125ml Syringe 1* J Tube (QFT)



P/N : 300108
Nemoto Press Duo & Press Pro &
Rempres Angio.

1*150ml Syringe 1* J Tube (QFT)





P/N : 400101
Bracco Empower CT & Empower CTA.

1*200ml Syringe 1*1500mm Coiled Tube 1* J Tube (QFT)



P/N : 400101S
Bracco Empower CT & Empower CTA.

1*200ml Syringe 1*1500mm Coiled Tube 1* Large Spike




P/N : 400101D
Bracco Empower CT& Empower CTA.

2*200ml Syringes 1* 1500mm Y Coiled Tube 2* J Tubes (QFTs)



P/N : 400101SD
Bracco Empower CT & Empower CTA.

2*200ml Syringes 1* 1500mm Y Coiled Tube 2*Spikes



P/N : 400102
Bracco Empower MR.

2*100ml Syringes 1* 2500mm T Coiled Tube
1* Small Spike 1*Large spike





P/N : 500101
Medtron Accutron CT.

1*200ml Syringe 1*1500mm Coiled Tube 1* J Tube (QFT)



P/N : 500101S
Medtron Accutron CT.

1*200ml Syringe 1*1500mm Coiled Tube 1* Large Spike




P/N : 500101D
Medtron Accutron CT-D.

2*200ml Syringes 1* 1500mm Y Coiled Tube 2* J Tubes (QFTs)




P/N : 500101SD
Medtron Accutron CT-D.

2*200ml Syringes 1* 1500mm Y Coiled Tube 2*Spikes




P/N : 500102
Medtron Accutron MR.

2*60ml Syringes 1* 2500mm T Coiled Tube
1* Small Spike 1*Large spike



P/N : 500103
Medtron Accutron HP.

1*200ml Syringe 1* J Tube (QFT)




P/N : 600101
Medrad Salient Single CT.

1*200ml Syringe 1*1500mm Coiled Tube 1* J Tube (QFT)




P/N : 600101S
Medrad Salient Single CT.

1*200ml Syringe 1*1500mm Coiled Tube 1* Large Spike




P/N : 600101D
Medrad Salient Dual CT.

2*200ml Syringes 1* 1500mm Y Coiled Tube 2* J Tubes (QFTs)



P/N : 600101SD
Medrad Salient Dual CT.

2*200ml Syringes 1* 1500mm Y Coiled Tube 2*Spikes



P/N : 600102
Medrad Avidia Angio.

1*150ml Syringes 1* J Tube (QFT)



P/N : 700101
Shenzhen ANKE ASA200 CT.

1*200ml Syringe 1*1500mm Coiled Tube 1* J Tube (QFT)



P/N : 700102
Shenzhen Sinopower-s CT.

1*200ml Syringe 1*1500mm Coiled Tube 1* J Tube (QFT)



P/N : 700102S
Shenzhen Sinopower-s CT.

1*200ml Syringe 1*1500mm Coiled Tube 1* Large Spike



P/N : 700102D
Shenzhen Sinopower-D CT.

2*200ml Syringes 1* 1500mm Y Coiled Tube 2* J Tubes (QFTs)



P/N : 700102SD
Shenzhen Sinopower-D CT.

2*200ml Syringes 1* 1500mm Y Coiled Tube 2*Spikes




P/N : 700103
Shenzhen Seacrown C20 CT.

1*200ml Syringe 1*1500mm Coiled Tube 1* J Tube (QFT)



P/N : 700103S
Shenzhen Seacrown C20 CT.

1*200ml Syringe 1*1500mm Coiled Tube 1* Large Spike



P/N : 700103D
Shenzhen Seacrown C22 CT.

2*200ml Syringes 1* 1500mm Y Coiled Tube 2* J Tubes (QFTs)



P/N : 700103SD
Shenzhen Seacrown C22 CT.

2*200ml Syringes 1* 1500mm Y Coiled Tube 2*Spikes



P/N : 700104
Shenzhen Seacrown C10 CT.

1*100ml Syringe 1*1500mm Coiled Tube 1* Large Spike



P/N : 700104S
Shenzhen Seacrown C10 CT.

1*100ml Syringe 1*1500mm Coiled Tube 1* Large Spike



P/N : 700104D
Shenzhen Seacrown C11 CT.

2*100ml Syringes 1* 1500mm Y Coiled Tube 2* J Tubes (QFTs)



P/N : 700104SD
Shenzhen Seacrown C11 CT.

2*100ml Syringes 1* 1500mm Y Coiled Tube 2*Spikes




P/N : 700105
Shenzhen SinoMRI AP/BP MR, Shenzhen Seacrown C60 MR.

2*65ml Syringes 1* 2500mm T Coiled Tube
1* Small Spike 1*Large spike



P/N : 700106
Shenzhen Sino Angio, Shenzhen Seacrown H15 Angio.

1*150ml Syringes 1* J Tube (QFT)



P/N : 800101
1500mm (60") Coiled Tube

1500mm Coiled Tube




P/N : 800102/800103
1500mm (60") / 2500mm (96") Y Coiled Tube With Valve.

1500mm Y Coiled Tube / 2500mm y Coiled Tube




P/N : 800104 / 800105 1500mm (60")/2500mm (96") T Coiled Tube With Valve

1500mm T Coiled Tube / 2500mm T Coiled Tube




P/N : 800106
All CT , ANGIO.

J Tube (QFT-Quick Fill Tube)




P/N : 800107/800108
Large Spike for Saline
Small Spike for Contrast agents

1*Large Spike 1* Small Spike



P/N : 800109
Used for Single Head Injector 12 hours & 24 hours injection.

1* Transfer Set with Single Chamber.



P/N : 800109D
Used for Dual Head Injector 12 hours & 24 hours injection.

2* Transfer set with single chamber 1*1500mm Y Coiled Tube



P/N : 800110
Used for transferring the medicine.
Clave connector.

1* Transfer Set with Clave Connector without chamber



P/N : SPT
Ulrich CT MRI Patient Tubing

SPT150	1500mm
SPT180	1800mm
SPT250	2500mm



P/N : HPT
1200psi High Pressure Tubing, Used for All Angio .



HPT050	500mm
HPT075	750mm
HPT100	1000mm
HPT120	1200mm
HPT150	1500mm
HPT180	1800mm
HPT200	2000mm
HPT250	2500mm



P/N : MPT
Ulrich CT MRI Daily use Pump Tubing.

1* Pump Tubing

Contacts

-  Service Hotline: 86-755-33557958
(9:30—12:00 12:00—18:00 Beijing Time)
-  Email: medis@medismed.com



MEDIS trademark belongs to Shenzhen Medis Medical Supplies Co., Ltd. Medis medical has the full copyright of Medis trademark, any other cannot use the trademark without the authorization, otherwise, medis medical will protect our legal right with the state law.

Version: 2022-V1

