



Product Name: Disposable High-pressure Syringes



## MRI Syringe Kits

### 65ml Syringe \*250cm T Connector \*Spikes

P/N : 100103

For use with Medrad MRI Injectors  
**Spectris MR Injectors**

#### Description:

- 2\*65ml Syringes
- 1\*2500mm(96") T Coiled Tube 350PSI
- 1\*Large spike 1\*Small Spike

#### Components / Materials

##### Luer Type

- Male Luer Connector
- Female Luer connector

##### Materials

- Polyethylene Terephthalate (PET)
- Synthetic Rubber
- Polycarbonate (PC)
- Polyvinyl Chloride (PVC)—TOTM
- / Polyvinyl Chloride (PVC)—DEHP
- Polyethylene (PE)

Exclusive Distributed by:

**Shenzhen Medis Medical Supplies Co., Ltd.**  
Rm1001, Block A, Building 3, Tianan Cloud Park, Bantian  
Street, Longgang District, Shenzhen 518129, China.  
Made in China  
[www.medismed.com](http://www.medismed.com)

#### Package / Sterilization

##### Package

- Blister package—Each kit
- Paper box package—Outer box
- Quantity per case—50kits
- Quantity per tray —1250kits
- Lid — Tyvek Medical paper
- Sterilization — EO GAS
- Gross Weight per case: 10kg
- Case size : 57\*39\*39cm

##### Shelf Life

3 Years

For more distributor information,  
please contact medismed direct  
[medis@medismed.com](mailto:medis@medismed.com)