



Product Name: Disposable High-pressure Syringes



STERILE EO



MRI Syringe Kits

65/115ml Syringe *250cm T Connector *Spikes

P/N : 100104

For use with Medrad MRI Injectors
Spectris Solaris MR Injectors

Description:

1*65ml Syringe 1*115ml Syringe
1*2500mm(96") T Coiled Tube 350PSI
1*Large spike 1*Small Spike

Components / Materials

Luer Type

Male Luer Connector
Female Luer connector

Materials

Polyethylene Terephthalate (PET)
Synthetic Rubber
Polycarbonate (PC)
Polyvinyl Chloride (PVC)—TOTM
/ Polyvinyl Chloride (PVC)—DEHP
Polyethylene (PE)

Exclusive Distributed by:

Shenzhen Medis Medical Supplies Co., Ltd.
Rm1001, Block A, Building 3, Tianan Cloud Park, Bantian
Street, Longgang District, Shenzhen 518129, China.
Made in China
www.medismed.com

Package / Sterilization

Package

Blister package—Each kit
Paper box package—Outer box
Quantity per case—50kits
Quantity per tray —1250kits
Lid — Tyvek Medical paper
Sterilization — EO GAS
Gross Weight per case: 11kg
Case size : 57*40*43cm

Shelf Life

3 Years

For more distributor information,
please contact medismed direct
medis@medismed.com

MEDIS
For Healthy Life