



Product Name: Disposable High-pressure Syringes



## Angio Syringe Kits 150ml Syringe \* J Tube (QFT)



P/N : 100108

For use with Medrad Angio Injectors  
**Mark VII & Avanta Injectors**

### Description:

1\*150ml Syringe

1\*J Tube (QFT)

### Components / Materials

Luer Type

Male Luer Connector

Female Luer connector

Materials

Polyethylene Terephthalate (PET)

Synthetic Rubber

Polycarbonate (PC)

Polyethylene (PE)

### Package / Sterilization

#### Package

Blister package—Each kit

Paper box package—Outer box

Quantity per case—50kits

Quantity per tray —1500kits

Lid — Tyvek Medical paper

Sterilization — EO GAS

Gross Weight per case: 6kg

Case size : 57\*31\*36cm

#### Shelf Life

3 Years

Exclusive Distributed by:

**Shenzhen Medis Medical Supplies Co., Ltd.**

Rm1001, Block A, Building 3, Tianan Cloud Park, Bantian Street, Longgang District, Shenzhen 518129, China.

Made in China

www.medismed.com

For more distributor information,  
please contact medismed direct

[medis@medismed.com](mailto:medis@medismed.com)

**MEDIS**  
For Healthy Life