



Product Name: Coiled Tubing



Braid Patient Tubing with valve
500mm--1800mm * 1200psi



P/N : HPT050—HPT180

For use with All series Angio Injectors

Description:

Length: 500mm --2500mm

Maximum Pressure: 1200psi

Inner diameter: 1.75mm---1.85mm

Outer diameter: 3.2mm ---3.3mm

Valve: With one valve

Latex : Free Latex

Components / Materials

Luer Type

Male Luer Connector

Female Luer connector

Materials

Polyvinyl Chloride (PVC)

Polyethylene (PE)

Acrylonitrile Butadine Styrene(ABS)

Package / Sterilization

Package

Paper Plastic package—Each kit

Paper box package—Outer box

Quantity per case—400kits

Lid — Tyvek Medical paper

Sterilization — EO GAS

Gross Weight per case: 13.5kg

Case size : 64*41*43cm

Shelf Life

3 Years

Exclusive Distributed by:

Shenzhen Medis Medical Supplies Co., Ltd.

Rm1001, Block A, Building 3, Tianan Cloud Park, Bantian Street, Longgang District, Shenzhen 518129, China.

Made in China

www.medismed.com

For more distributor information,
please contact medismed direct

medis@medismed.com

MEDIS
For Healthy Life