



Product Name: Disposable High-pressure Syringes



## Dual CT Syringe Kits

### 200ml Syringe \*150cm Y Connector \*J Tube

P/N : 200101D

For use with LF CT Injectors


**CT9000 & CT9000 ADV, Optivantage DH, OptiOne CT**

Description:

2\*200ml Syringes

1\*1500mm(60") Y Coiled Tube 350PSI

2\*J Tubes



### Components / Materials

Luer Type

Male Luer Connector

Female Luer connector

Materials

Polypropylene(PP)

Synthetic Rubber

Polycarbonate (PC)

Polyvinyl Chloride (PVC)—TOTM

/ Polyvinyl Chloride (PVC)—DEHP

Polyethylene (PE)

Exclusive Distributed by:

**Shenzhen Medis Medical Supplies Co., Ltd.**

Rm1001, Block A, Building 3, Tianan Cloud Park, Bantian Street, Longgang District, Shenzhen 518129, China.

Made in China

[www.medismed.com](http://www.medismed.com)

### Package / Sterilization

#### Package

Blister package—Each kit

Paper box package—Outer box

Quantity per case—20kits

Quantity per tray —500kits

Lid — Tyvek Medical paper

Sterilization — EO GAS

Gross Weight per case: 7.5kg

Case size : 49\*36\*38cm

#### Shelf Life

3 Years

For more distributor information, please contact medismed direct

[medis@medismed.com](mailto:medis@medismed.com)