




Product Name: Disposable High-pressure Syringes




## Dual CT Syringe Kits 200ml Syringe \*150cm Y Connector \*Spikes

P/N : 200101SD

For use with LF CT Injectors  
**CT9000 & CT9000 ADV, Optivantage DH, OptiOne CT**

### Description:

- 2\*200ml Syringes
- 1\*1500mm(60") Y Coiled Tube 350PSI
- 1\* Large Spike 1\* Small Spike



### Components / Materials

#### Luer Type

- Male Luer Connector
- Female Luer connector

#### Materials

- Polypropylene(PP)
- Synthetic Rubber
- Polycarbonate (PC)
- Polyvinyl Chloride (PVC)—TOTM
- / Polyvinyl Chloride (PVC)—DEHP
- Polyethylene (PE)

Exclusive Distributed by:

**Shenzhen Medis Medical Supplies Co., Ltd.**  
Rm1001, Block A, Building 3, Tianan Cloud Park, Bantian  
Street, Longgang District, Shenzhen 518129, China.  
Made in China  
[www.medismed.com](http://www.medismed.com)

### Package / Sterilization

#### Package

- Blister package—Each kit
- Paper box package—Outer box
- Quantity per case—20kits
- Quantity per tray —500kits
- Lid — Tyvek Medical paper
- Sterilization — EO GAS
- Gross Weight per case: 7.5kg
- Case size : 49\*36\*38cm

#### Shelf Life

3 Years

For more distributor information,  
please contact medismed direct  
[medis@medismed.com](mailto:medis@medismed.com)