



Product Name: Disposable High-pressure Syringes



STERILE EO



## MRI Syringe Kits

60ml Syringe \*250cm T Connector  
\*Spikes

P/N : 200102

For use with LF MR Injectors  
**Optistar LE, Elite MR**

### Description:

2\*60ml Syringes  
1\*2500mm(96") T Coiled Tube 350PSI  
1\*Large Spike 1\* Small Spike

### Components / Materials

#### Luer Type

Male Luer Connector  
Female Luer connector

#### Materials

Polyethylene Terephthalate (PET)  
Synthetic Rubber  
Polycarbonate (PC)  
Polyvinyl Chloride (PVC)—TOTM  
/ Polyvinyl Chloride (PVC)—DEHP  
Polyethylene (PE)

Exclusive Distributed by:

**Shenzhen Medis Medical Supplies Co., Ltd.**  
Rm1001, Block A, Building 3, Tianan Cloud Park, Bantian  
Street, Longgang District, Shenzhen 518129, China.  
Made in China  
[www.medismed.com](http://www.medismed.com)

### Package / Sterilization

#### Package

Blister package—Each kit  
Paper box package—Outer box  
Quantity per case—50kits  
Quantity per tray —1250kits  
Lid — Tyvek Medical paper  
Sterilization — EO GAS  
Gross Weight per case: 10kg  
Case size : 57\*39\*39cm

#### Shelf Life

3 Years

For more distributor information,  
please contact medismed direct  
[medis@medismed.com](mailto:medis@medismed.com)

**MEDIS**  
For Healthy Life