



Product Name: Disposable High-pressure Syringes



Angio Syringe Kits 150ml Syringe* J Tube (QFT)



STERILE EO



P/N : 200103

For use with LF Angio Injectors
Angiomat Illumena Neo Angio

Description:

1*150ml Syringe

1*J Tube

Components / Materials

Luer Type

Male Luer Connector

Female Luer connector

Materials

Polyethylene Terephthalate (PET)

Synthetic Rubber

Polycarbonate (PC)

Polyethylene (PE)

Package / Sterilization

Package

Blister package—Each kit

Paper box package—Outer box

Quantity per case—50kits

Quantity per tray —1500kits

Lid — Tyvek Medical paper

Sterilization — EO GAS

Gross Weight per case: 7kg

Case size : 57*31*36cm

Shelf Life

3 Years

Exclusive Distributed by:

Shenzhen Medis Medical Supplies Co., Ltd.

Rm1001, Block A, Building 3, Tianan Cloud Park, Bantian Street, Longgang District, Shenzhen 518129, China.

Made in China

www.medismed.com

For more distributor information,
please contact medismed direct

medis@medismed.com

MEDIS
For Healthy Life