



Product Name: Disposable High-pressure Syringes



Single CT Syringe Kits 200ml Syringe *150cm Connector *Spike



STERILE EO



P/N : 300102S

For use with Nemoto CT Injectors
A25 / A60 / A300 CT

Description:

- 1*200ml Syringe
- 1*1500mm(60") Coiled Tube 350PSI
- 1*Large Spike



Components / Materials

Luer Type

- Male Luer Connector
- Female Luer connector

Materials

- Polyethylene Terephthalate (PET)
- Synthetic Rubber
- Polycarbonate (PC)
- Polyvinyl Chloride (PVC)—TOTM
- / Polyvinyl Chloride (PVC)—DEHP
- Polyethylene (PE)

Exclusive Distributed by:

Shenzhen Medis Medical Supplies Co., Ltd.

Rm1001, Block A, Building 3, Tianan Cloud Park, Bantian Street, Longgang District, Shenzhen 518129, China.

Made in China

www.medismed.com

Package / Sterilization

Package

- Blister package—Each kit
- Paper box package—Outer box
- Quantity per case—50kits
- Quantity per tray —1250kits
- Lid — Tyvek Medical paper
- Sterilization — EO GAS
- Gross Weight per case: 9kg
- Case size : 57*31*37cm

Shelf Life

3 Years

For more distributor information,
please contact medismed direct

medis@medismed.com

MEDIS
For Healthy Life