




Product Name: Disposable High-pressure Syringes



## Dual CT Syringe Kits 200ml/100ml Syringe \*150cm Y Connector\*J Tubes

P/N : 300103



For use with Nemoto CT Dual shot Injectors  
**GX-V & Alpha 7 & B100 & B200 CT**

### Description:

1\*200ml Syringe 1\*100ml Syringe  
1\*1500mm(60") Y Coiled Tube 350PSI  
2\*J Tubes

### Components / Materials

#### Luer Type

Male Luer Connector  
Female Luer connector

#### Materials

Polyethylene Terephthalate (PET)  
Synthetic Rubber  
Polycarbonate (PC)  
Polyvinyl Chloride (PVC)—TOTM  
/ Polyvinyl Chloride (PVC)—DEHP  
Polyethylene (PE)

Exclusive Distributed by:

**Shenzhen Medis Medical Supplies Co., Ltd.**  
Rm1001, Block A, Building 3, Tianan Cloud Park, Bantian  
Street, Longgang District, Shenzhen 518129, China.  
Made in China  
[www.medismed.com](http://www.medismed.com)

### Package / Sterilization

#### Package

Blister package—Each kit  
Paper box package—Outer box  
Quantity per case—20kits  
Quantity per tray —500kits  
Lid — Tyvek Medical paper  
Sterilization — EO GAS  
Gross Weight per case: 9kg  
Case size : 49\*32\*33cm

#### Shelf Life

3 Years

For more distributor information,  
please contact [medismed direct](mailto:medismed@medismed.com)  
[medismed@medismed.com](mailto:medismed@medismed.com)