



Product Name: Disposable High-pressure Syringes



Dual CT Syringe Kits

100ml/60ml Syringe


*150cm Y Connector*J Tubes

P/N : 300105

For use with Nemoto CT Dual shot Injectors
Alpha 7 CT

Description:

1*100ml Syringe 1*60ml Syringe
1*1500mm(60") Y Coiled Tube 350PSI
2*J Tubes



Components / Materials

Luer Type

Male Luer Connector
Female Luer connector

Materials

Polyethylene Terephthalate (PET)
Synthetic Rubber
Polycarbonate (PC)
Polyvinyl Chloride (PVC)—TOTM
/ Polyvinyl Chloride (PVC)—DEHP
Polyethylene (PE)

Exclusive Distributed by:

Shenzhen Medis Medical Supplies Co., Ltd.
Rm1001, Block A, Building 3, Tianan Cloud Park, Bantian
Street, Longgang District, Shenzhen 518129, China.
Made in China
www.medismed.com

Package / Sterilization

Package

Blister package—Each kit
Paper box package—Outer box
Quantity per case—50kits
Quantity per tray —1250kits
Lid — Tyvek Medical paper
Sterilization — EO GAS
Gross Weight per case: 8kg
Case size : 49*32*33cm

Shelf Life

3 Years

For more distributor information,
please contact [medismed direct](mailto:medismed@medismed.com)
medismed@medismed.com