



Product Name: Disposable High-pressure Syringes



STERILE EO



## Angio Syringe Kits 125ml Syringe\*J Tube

P/N : 300107

For use with Nemoto Angio Injectors  
**120S & Press Pro Angio injectors**

### Description:

1\*125ml Syringe

1\*J Tube

### Components / Materials

Luer Type

Male Luer Connector

Female Luer connector

Materials

Polypropylene(PP)

Synthetic Rubber

Polycarbonate (PC)

Polyethylene (PE)

### Package / Sterilization

#### Package

Blister package—Each kit

Paper box package—Outer box

Quantity per case—50kits

Quantity per tray —1250kits

Lid — Tyvek Medical paper

Sterilization — EO GAS

Gross Weight per case: 6kg

Case size : 57\*31\*36cm

#### Shelf Life

3 Years

Exclusive Distributed by:

**Shenzhen Medis Medical Supplies Co., Ltd.**

Rm1001, Block A, Building 3, Tianan Cloud Park, Bantian Street, Longgang District, Shenzhen 518129, China.

Made in China

[www.medismed.com](http://www.medismed.com)

For more distributor information,  
please contact medismed direct  
[medis@medismed.com](mailto:medis@medismed.com)

**MEDIS**  
For Healthy Life