




Product Name: Disposable High-pressure Syringes




Dual CT Syringe Kits 200ml Syringes *150cm Y Connector *Spikes

P/N : 400101SD

For use with Bracco CT Injectors
EZEM Empower CT & CTA-D

Description:

- 2*200ml Syringes
- 1*1500mm(60") Y Coiled Tube 350PSI
- 2*Spikes



Components / Materials

Luer Type

- Male Luer Connector
- Female Luer connector

Materials

- Polyethylene Terephthalate (PET)
- Synthetic Rubber
- Polycarbonate (PC)
- Polyvinyl Chloride (PVC)—TOTM
- / Polyvinyl Chloride (PVC)—DEHP
- Polyethylene (PE)

Exclusive Distributed by:

Shenzhen Medis Medical Supplies Co., Ltd.

Rm1001, Block A, Building 3, Tianan Cloud Park, Bantian Street, Longgang District, Shenzhen 518129, China.

Made in China

www.medismed.com

Package / Sterilization

Package

- Blister package—Each kit
- Paper box package—Outer box
- Quantity per case—20kits
- Quantity per tray —500kits
- Lid — Tyvek Medical paper
- Sterilization — EO GAS
- Gross Weight per case: 8kg
- Case size : 49*36*38cm

Shelf Life

3 Years

For more distributor information,
please contact [medismed direct](mailto:medismed@medismed.com)
medis@medismed.com

MEDIS
For Healthy Life