



Product Name: Disposable High-pressure Syringes



STERILE EO



Single CT Syringe Kits

190ml Syringe *150cm Connector *J Tube

P/N : 600101

For use with Medrad CT Injectors
Imaxeon Salient

Description:

- 1*190ml Syringe
- 1*1500mm(60") Coiled Tube 350PSI
- 1*J Tube



Components / Materials

Luer Type

- Male Luer Connector
- Female Luer connector

Materials

- Polyethylene Terephthalate (PET)
- Synthetic Rubber
- Polycarbonate (PC)
- Polyvinyl Chloride (PVC)—TOTM
- / Polyvinyl Chloride (PVC)—DEHP
- Polyethylene (PE)

Exclusive Distributed by:

Shenzhen Medis Medical Supplies Co., Ltd.
Rm1001, Block A, Building 3, Tianan Cloud Park, Bantian
Street, Longgang District, Shenzhen 518129, China.
Made in China
www.medismed.com

Package / Sterilization

Package

- Blister package—Each kit
- Paper box package—Outer box
- Quantity per case—50kits
- Quantity per tray —1250kits
- Lid — Tyvek Medical paper
- Sterilization — EO GAS
- Gross Weight per case: 10kg
- Case size : 57*39*39cm

Shelf Life

3 Years

For more distributor information,
please contact medismed direct
medis@medismed.com

MEDIS
For Healthy Life