



Product Name: Disposable High-pressure Syringes



## Angio Syringe Kits 150ml Syringe\*J Tube



STERILE EO



P/N : 600102

For use with Medrad Angio Injectors  
**Imaxeon Avidia Angio**

### Description:

1\*150ml Syringe

1\*J Tube

### Components / Materials

Luer Type

Male Luer Connector

Female Luer connector

Materials

Polyethylene Terephthalate (PET)

Synthetic Rubber

Polycarbonate (PC)

Polyvinyl Chloride (PVC)—TOTM

/ Polyvinyl Chloride (PVC)—DEHP

Polyethylene (PE)

Exclusive Distributed by:

**Shenzhen Medis Medical Supplies Co., Ltd.**

Rm1001, Block A, Building 3, Tianan Cloud Park, Bantian Street, Longgang District, Shenzhen 518129, China.

Made in China

[www.medismed.com](http://www.medismed.com)

### Package / Sterilization

#### Package

Blister package—Each kit

Paper box package—Outer box

Quantity per case—50kits

Quantity per tray —1250kits

Lid — Tyvek Medical paper

Sterilization — EO GAS

Gross Weight per case: 6kg

Case size : 57\*31\*36cm

#### Shelf Life

3 Years

For more distributor information,  
please contact [medis@medismed.com](mailto:medis@medismed.com)

[medis@medismed.com](mailto:medis@medismed.com)