



Product Name: Disposable High-pressure Syringes



## Dual CT Syringe Kits


### 200ml Syringes \*150cm Y Connector \*J Tubes

P/N : 700102D

For use with SINO CT Injectors  
**Sinopower-D CT**

#### Description:

2\*200ml Syringes  
1\*1500mm(60") Y Coiled Tube 350PSI  
2\*J Tubes



#### Components / Materials

Luer Type

Male Luer Connector  
Female Luer connector

Materials

Polyethylene Terephthalate (PET)  
Synthetic Rubber  
Polycarbonate (PC)  
Polyvinyl Chloride (PVC)—TOTM  
/ Polyvinyl Chloride (PVC)—DEHP  
Polyethylene (PE)

Exclusive Distributed by:

**Shenzhen Medis Medical Supplies Co., Ltd.**  
Rm1001, Block A, Building 3, Tianan Cloud Park, Bantian  
Street, Longgang District, Shenzhen 518129, China.  
Made in China  
[www.medismed.com](http://www.medismed.com)

#### Package / Sterilization

##### Package

Blister package—Each kit  
Paper box package—Outer box  
Quantity per case—50kits  
Quantity per tray —1250kits  
Lid — Tyvek Medical paper  
Sterilization — EO GAS  
Gross Weight per case: 8kg  
Case size : 49\*36\*38cm

##### Shelf Life

3 Years

For more distributor information,  
please contact medismed direct  
[medis@medismed.com](mailto:medis@medismed.com)