



Product Name: Disposable High-pressure Syringes



STERILE EO



## Single CT Syringe Kits 100ml Syringe \*150cm Connector \*J Tube

P/N : 700104

For use with Seacrown CT Injectors  
**Zenith C10 CT**

### Description:

- 1\*200ml Syringe
- 1\*1500mm(60") Coiled Tube 350PSI
- 1\*J Tube

### Components / Materials

#### Luer Type

- Male Luer Connector
- Female Luer connector

#### Materials

- Polyethylene Terephthalate (PET)
- Synthetic Rubber
- Polycarbonate (PC)
- Polyvinyl Chloride (PVC)—TOTM
- / Polyvinyl Chloride (PVC)—DEHP
- Polyethylene (PE)

Exclusive Distributed by:

**Shenzhen Medis Medical Supplies Co., Ltd.**  
Rm1001, Block A, Building 3, Tianan Cloud Park, Bantian  
Street, Longgang District, Shenzhen 518129, China.  
Made in China  
[www.medismed.com](http://www.medismed.com)

### Package / Sterilization

#### Package

- Blister package—Each kit
- Paper box package—Outer box
- Quantity per case—50kits
- Quantity per tray —1250kits
- Lid — Tyvek Medical paper
- Sterilization — EO GAS
- Gross Weight per case: 9kg
- Case size : 57\*39\*39cm

#### Shelf Life

3 Years

For more distributor information,  
please contact medismed direct  
[medis@medismed.com](mailto:medis@medismed.com)