



Product Name: Disposable High-pressure Syringes



STERILE EO



Dual CT Syringe Kits

100ml Syringes *150cm Y Connector *Spikes

P/N : 700104SD

For use with Seacrown CT Injectors

Zenith C11

Description:

2*100ml Syringes

1*1500mm(60") Y Coiled Tube 350PSI

2*Spikes

Components / Materials

Luer Type

Male Luer Connector

Female Luer connector

Materials

Polyethylene Terephthalate (PET)

Synthetic Rubber

Polycarbonate (PC)

Polyvinyl Chloride (PVC)—TOTM

/ Polyvinyl Chloride (PVC)—DEHP

Polyethylene (PE)

Exclusive Distributed by:

Shenzhen Medis Medical Supplies Co., Ltd.

Rm1001, Block A, Building 3, Tianan Cloud Park, Bantian Street, Longgang District, Shenzhen 518129, China.

Made in China

www.medismed.com

Package / Sterilization

Package

Blister package—Each kit

Paper box package—Outer box

Quantity per case—20kits

Quantity per tray —500kits

Lid — Tyvek Medical paper

Sterilization — EO GAS

Gross Weight per case: 8kg

Case size : 49*36*38cm

Shelf Life

3 Years

For more distributor information, please contact medismed direct

medis@medismed.com