

CT

PRODUCT DATA SHEET



MRI

Product Name: Coiled Tubing



**Y connector with check valve
1500mm/2500mm * 350psi**



STERILE EO



P/N : 800102 / 800103

For use with All series CT MRI Injectors

Description:

Length: 1500mm (2500mm)

Maximum Pressure: 350psi

Inner diameter: 1.75mm---1.85mm

Outer diameter: 3.2mm ---3.3mm

Valve: With one valve

Latex : Free Latex

Components / Materials

Luer Type

Male Luer Connector

Female Luer connector

Materials

Polyvinyl Chloride (PVC)

Polyethylene (PE)

Acrylonitrile Butadiene Styrene (ABS)

Package / Sterilization

Package

Paper Plastic package—Each kit

Paper box package—Outer box

Quantity per case—400kits

Lid — Tyvek Medical paper

Sterilization — EO GAS

Gross Weight per case: 13.5kg

Case size : 64*41*43cm

Shelf Life

3 Years

Exclusive Distributed by:

Shenzhen Medis Medical Supplies Co., Ltd.

Rm1001, Block A, Building 3, Tianan Cloud Park, Bantian Street, Longgang District, Shenzhen 518129, China.

Made in China

www.medismed.com

For more distributor information,
please contact medismed direct

medis@medismed.com

MEDIS
For Healthy Life